

Heat activation of TRPV3 on Nanion's Patchliner

Nanion Technologies GmbH, Munich

Summary

The transient receptor potential cation channel, subfamily V, member 3 (TRPV3), is a ligand-gated, non-selective cation channel first described in 2002^{1,2,3}. It exhibits 43% sequence identity to TRPV1¹. Although TRPV3 has been detected immunologically in the CNS and suggested to be often co-localized with TRPV1³, it is found more robustly in keratinocytes in skin and, given its threshold for temperature activation of $>34^{\circ}\text{C}^{1,2,3}$, it has been speculated that TRPV3 may act in co-operation with sensory afferents to perceive warmth and signal elevated temperature. TRPV3 can also be activated by the ligand 2-Aminoethoxydiphenyl borate (2-APB)⁴. The TRPV3 channel, along with other TRPV channels, may play an important role in chronic pain and, therefore, is receiving increasing attention as a potential therapeutic candidate for the treatment of chronic pain.

Here we present data collected on a 4-channel Patchliner with temperature control showing the potential use of the Patchliner to record TRPV3 currents activated by 2-APB or heat. As previously reported, TRPV3 currents sensitize to repeated applications of 2-APB⁴ or heat¹, a phenomenon we also observed. At low concentrations of 2-APB, the currents were primarily outwardly rectifying but at higher concentrations and with prolonged exposure they often became dual rectifying (data not shown). This is also in good agreement with the literature⁴. In contrast, the temperature-activated responses were always outwardly rectifying with little inward current. The inward currents activated by 2-APB could be blocked by ruthenium red (RR) as expected⁴.

Results

Current responses of an individual cell to 200 ms voltage ramps (-100 mV to $+100\text{ mV}$) and activation by application of increasing concentrations of 2-APB are shown in Figure 1. A concentration response curve (using amplitude at $+90\text{ mV}$) revealed an EC_{50} for 2-APB activation of $61.1 \pm 7.6\text{ }\mu\text{M}$ ($n = 11$), in good agreement with the literature^{4,5}. Interestingly, as previously reported^{1,4}, TRPV3 currents displayed sensitization to repeated application of 2-APB or heat (Figure 1).

Figure 2 shows the activation of TRPV3 channels by increasing temperature. Outwardly rectifying currents started to activate at 38°C and increased in amplitude as temperatures increased, up to 54°C , in good agreement with the literature^{1,2,3}. To rule out the involvement of TRPV1 a control was performed using capsaicin. Currents activated by $100\text{ }\mu\text{M}$ 2-APB were not activated by $1\text{ }\mu\text{M}$ capsaicin (Figure 2B) as previously reported¹.

A full concentration response curve (CRC) to RR at 0.1 , 1 , 3 and $10\text{ }\mu\text{M}$ was performed (Figure 2C,D) using $100\text{ }\mu\text{M}$ 2-APB as the




Figure 1: **A** Activation of TRPV3 by increasing concentrations of 2-APB. **B** Concentration response curve for 2-APB activation, $\text{EC}_{50} = 61.1 \pm 7.6\text{ }\mu\text{M}$ ($n = 11$). **C & D** Repeated exposure of TRPV3 currents to 2-APB (**C**) or temperature (**D**) caused sensitization of currents (lines indicate 2-APB or heat exposure).




Figure 2: **A** TRPV3 was activated by increasing temperature. **B** TRPV3 was activated by 2-APB but not capsaicin. **C** The inward TRPV3 current was blocked by RR in a concentration-dependent manner but not the outward current (not shown). **D** CRC of the inward current (at -50 mV) revealed an $IC_{50} = 2.4 \pm 0.4 \mu M$ ($n = 4$).

activator. As expected⁴, RR blocked the inward current with an IC_{50} of $2.4 \pm 0.4 \mu M$ ($n = 4$). The outward current was not blocked by RR (not shown).

In summary, TRPV3 receptors stably expressed in HEK293 cells can be reliably activated by moderate ($\geq 38^\circ C$) to hot temperatures (up to $54^\circ C$), or by 2-APB in a concentration-dependent manner. The data shown here agrees well with published literature using conventional patch clamp electrophysiology to study TRPV3¹⁻⁵. Therefore, the Patchliner provides a viable, higher throughput alternative to conventional patch clamp for the discovery of active TRPV3 lead compounds with a suitable drug profile. What is more, using the unique feature of the Patchliner which allows rapid and transient activation of TRP channels by heating the solution inside the pipette, compounds that block the ligand response could be separated from those that block the heat-activated response.

References

1. Peier, A.M., *et al.*, 2002. Science. 296: 2046-2049.
2. Xu, H., *et al.*, 2002. Nature. 418: 181-186.
3. Smith, G.D., *et al.*, 2002. Nature. 418: 186-190.
4. Chung, M.K., *et al.*, 2004. J. Neurosci. 24: 5177-5182.
5. Hu *et al.*, 2009. PNAS. 106: 1626-1631.

Methods

Cells

Precislon hTRPV3-HEK recombinant cell line (CYL3065) was kindly provided by EMD Millipore, USA.

Electrophysiology measurements

Cells were cultured and harvested according to Nanion's

standard cell culture protocol. Cells were resuspended in external recording solution and stored in the CellHotel of the Patchliner before being dispensed into each well of the NPC-16 chip. Internal and external solution compositions are available upon request. Whole cell patch clamp recordings were conducted according to Nanion's standard procedure for the Patchliner. Currents were elicited by 200 ms voltage ramps from -100 mV to +100 mV, $V_{hold} = -20$ mV. For heat activation of channels, 100 μl of external solution was heated in the pipette until desired temperature was reached (within 30 s) and rapidly applied to the cell at a speed of 10 $\mu l/s$ (N.B. a control application of 50 μM 2-APB at room temperature was performed at the start to ensure expression of TRPV3 in each cell prior to heat application). 2-APB was diluted in external solution at the indicated concentrations and applied at room temperature for approximately 35 s before wash with external solution. RR was diluted in external solution at the indicated concentrations and pre-incubated for at least 30 s before co-application with 100 μM 2-APB.

